

# Muse iQ CROS System

---

*Hearing solutions for  
single-sided hearing loss*



Call today to learn more  
877-826-2941





# Made for single-sided hearing loss

---

With our new line of **Muse™ iQ** hearing aids, Starkey® has a solution for those with single-sided hearing loss.

Starkey can help you surround yourself with the power of Starkey Sound.

*Starkey Sound™ uses cutting-edge technology to deliver true listening enjoyment.*



Muse™ iQ



# The Muse iQ CROS System

---



*The new Muse iQ CROS System was developed to give more people solutions for better hearing.*

---

Now, if you have single-sided hearing loss you'll be able to hear the world around you in full, pristine sound.



**Muse iQ CROS  
System benefits:**

- 1.** Allows you to hear from both sides — with devices so small and discreet they're nearly invisible
- 2.** Improves listening in noise
- 3.** Reduces background noise

# Customized solutions



## CROS solution

If you're unable to hear in one ear and have normal hearing in the other ear, you may benefit from a CROS solution. Finally, you'll be able to hear sounds in your weaker-hearing ear by way of your healthy ear.

### How it works

To help you better engage in all types of listening environments, our CROS solution has a microphone that picks up sounds and voices from your weaker-hearing ear, then transmits them to a hearing aid receiver that's fitted on your good ear. No matter your position at the dinner table, or if sounds originate from your weaker-hearing side, you won't miss important sounds around you.



**Also available with a**



## BiCROS solution

A BiCROS solution is also available and can benefit those with little to no hearing in one of their ears, and a hearing loss in their better ear.

### How it works

This option sends sound from a microphone, placed in your ear with little to no hearing, to a hearing aid receiver placed in your better-hearing ear. This allows you to hear from both sides. Your hearing care professional can help you determine what solution will best fit your needs and lifestyle.

# Made for *convenient connectivity*

With more and more wireless products on the market, making connections between devices can be challenging. To navigate listening environments with high wireless traffic, our 900sync™ technology provides a powerful solution. Its high bit-rate streaming allows for ear-to-ear phone streaming and wireless capabilities for single-sided hearing loss.



# SURFLINK®

*Our Starkey wireless hearing aids wouldn't be complete without our popular SurfLink accessories.*



Our compact cell phone transmitter, **SurfLink® Mobile 2**, is an assistive listening device, media streamer and hearing aid remote all rolled into one.



**SurfLink Mini Mobile** is a small and simple cell phone streaming solution. With this device, you can connect any Bluetooth®-ready device to your Muse iQ hearing aids.



**SurfLink Media 2** is a set-and-forget media streaming solution that connects to TVs, MP3 players and more, to wirelessly stream audio directly to your hearing aids.



Our **SurfLink Remote** can be used to adjust volume, change memory modes, go in and out of streaming mode and much more — all with the push of a button.



Our discreet, lightweight **SurfLink Remote Microphone 2** is worn by a conversation partner to aid patients in one-on-one conversations in difficult listening environments.



# Made to deliver the very best

You want

*To hear comfortably in  
challenging environments*

---

*To easily adapt to  
new sounds*

---

*Better wireless  
streaming*

---

*A comfortable, personalized  
listening experience*

---

*No buzzing or whistling*

---

*Durable, dependable  
hearing aids*

---

*Music the way it was  
meant to be heard*

---

*Customizable tinnitus relief*

---

Features may vary by technology level.



## Muse iQ delivers

NEW

**Acuity™ Immersion Directionality**, designed to restore natural benefits provided by the ear and improve speech audibility in difficult listening situations.

This feature speeds your transition to hearing aids by gradually adjusting settings to allow your brain to get used to new sounds.

Technology that provides consistent wireless performance while streaming cell phones, TV, music and other media with **SurfLink Mobile 2**, **SurfLink Mini Mobile** and other **SurfLink** accessories.

Capabilities that speed your transition to hearing aids with gradual adjustments, distortion-free listening in noise, and ultimate clarity for soft sounds.

Starkey's best-in-class feedback cancellation system providing feedback-free, comfortable listening all-day long.

**Surface™ NanoShield**, our pioneering water, wax and moisture repellent system, to protect and ensure durability and dependability.

A design that allows you to experience music in a whole new way. Now you can hear every note the way the artist intended — with pure, refined sound quality.

Advanced **Multiflex Tinnitus Technology** integrated into Muse iQ to bring relief to those who suffer from ringing in the ears.

# Made for everywhere you are


***From the beginning,  
our hearing solutions  
have enabled people to  
hear better in challenging  
environments and  
live a better life.***

You want to hear all the subtleties of life: the inflection in the voice of a spouse or the nuances in the notes of a favorite song. Whether it's one-on-one or in the most crowded places, you'll hear the world around you.

What you need will depend on how busy and active you are, and what types of environments you frequent.



# Listening environments



Technology Levels

*What is your  
hearing  
lifestyle?*

Our lineup of products enhances many different listening environments. Talk to your hearing professional about which options fit your lifestyle best.

*Compare the listening  
environments on the  
chart inside.*

Made for  
your style

*Styles and colors*



Muse iQ

Pick your favorite color among our new RIC and BTE styles, both with tinnitus and wireless CROS and BiCROS solutions.

Consult your hearing care professional and start experiencing Muse iQ's discreet, sleek styles for yourself.

## COLOR GUIDE



Black



Slate



Sterling



Bright White  
with Sterling



Espresso



Bronze



Champagne



*Muse iQ is available  
in a variety of styles  
and colors for your  
ideal fit.*



# Hear the world in surround sound.

---

***Starkey delivers the best possible sound to your better ear so you don't have to worry about repositioning yourself for conversations or missing important sounds.***

Single-sided hearing loss doesn't have to get in the way of your lifestyle or your activities. With the Muse iQ CROS System, you can experience sound from all sides — never missing a single moment that brings you joy. Because when you hear better you live better.





Muse iQ



Made so  
the world  
may hear



***Our Founder and CEO,  
William F. Austin, has always  
held a strong belief that the  
gift of better hearing results  
in a fuller and richer life.***

We prove our dedication to this philosophy daily, through our support of Starkey Hearing Foundation. And for every Starkey hearing aid purchased, we help someone in need.

The gift of hearing opens up new worlds of opportunity. It connects individuals to life and helps them realize they can accomplish more than they ever thought possible. Starkey Hearing Foundation uses hearing as a vehicle to change lives around the world.



*Made to help millions.*



**100+**  
COUNTRIES



**1.1** MILLION PEOPLE  
IN NEED SERVED



**COUNTLESS**  
LIVES CHANGED

Spread hope from ear to ear. To learn more or give the gift of hearing, visit [StarkeyHearingFoundation.org](http://StarkeyHearingFoundation.org).

# Questions?

Call today to speak with  
a hearing care expert

877-826-2941



[www.Starkey.com](http://www.Starkey.com)

Global Headquarters  
6700 Washington Ave. S.  
Eden Prairie, MN 55344  
1.800.328.8602